

DOCUMENTATION LABEL AND GUN



PRODUCT FEATURES

ZEISS labels are produced under ISO 13485 systems. Labels are produced for use, Ethylene Steam and Hydrogen Labels can be peeled off after it can be kept or department records for proof of process. Labels Type1 ACCORDING TO EN ISO11140-1 STANDARD

CODE	PRODUCT NAME	FEATURE
ZD 10	EO Documation Label	600 pieces
ZD 20	Steam Documation Label	600 pieces
ZD 30	Hidrogen Peroxide Documation Label	600 pieces
ZDT1	Documation Label Gun	1