

ETHYLENE OXIDE GAS CARTRIDGE/AMPOULE KIT



- EO GAS CARTRIDGE
AMPOULE KIT**
(The contents of each Kit)
- EO Gas Cartridge/Ampoule
 - Sterilization Bag
 - Humidity Chip
 - Dosimeters
 - Clip

They are the sterilization kits, operated with the gas diffusion technology which sterilize all kinds of plastic, rubber, sensitive, metal, mechanical, medical and surgical materials whether they are resistant to heat or not.



AMPOULE KIT					
MODEL	ACTIVE MATERIAL ETHYLENE OXIDE	WRAPPED	STERILIZATION BAG SIZE	NUMBER OF KIT IN ONE BOX	NET CONTENT
ZQ-11	100 %	Poliüretan	60x90	100	11 gr
ZQ-15	100 %	Poliüretan	60x90	100	15 gr
ZQ-18	100 %	Poliüretan	60x90	100	18 gr
ZQ-20	100 %	Poliüretan	60x90	100	20 gr
ZQ-22	100 %	Poliüretan	60x90	100	22 gr
ZQ-24	100 %	Poliüretan	94x135	100	24 gr
ZQ-30	100 %	Poliüretan	94x135	100	30 gr

CARTRIDGE KIT					
MODEL	ACTIVE MATERIAL ETHYLENE OXIDE	WRAPPED	STERILIZATION BAG SIZE	NUMBER OF KIT IN ONE BOX	NET CONTENT
ZS-11	100 %	Poliüretan	60x90	50	11 gr
ZS-15	100 %	Poliüretan	60x90	50	15 gr
ZS-18	100 %	Poliüretan	60x90	50	18 gr
ZS-20	100 %	Poliüretan	60x90	50	20 gr
ZS-22	100 %	Poliüretan	60x90	50	22 gr