

# ETHYLENE OXIDE ALUMINUM GAS CARTRIDGES



**%100 EO**

Med Tip provides ZEISS® brand ethylene oxide in aluminum cartridges with controlled quantities and packed under strict safety regulation that meet the requirements of the variety of existing equipment in the market. This system provides safety in the sterilization areas, as the cartridge can only be punctured inside the sterilizer for which it is designed.

## Features

- Ethylene oxide of high pureness (100% ethylene oxide).
- Single-use Cartridges.
- High effectiveness as a bactericide, fungicide and antivirus.
- Compatible size with any sterilizer in the market.

• **3M Sterivac® - Amsco® - Hogner® - Del Giudice® Biomedica® - and many others.**

- Do not contain CFC (freon), do not affect the ozone layer.

CODE	MEASURE		CAPACITY	
	Diameter	Height	Grams	Oz.
ZL-100	38 mm	165 mm	100 gr	3,52 oz
ZL-127	38 mm	165 mm	127 gr	4,47 oz
ZL-130	50 mm	133 mm	130 gr	4,60 oz
ZL-134	38 mm	165 mm	134 gr	4,72 oz
ZL-170	45 mm	178 mm	170 gr	5,99 oz
ZL-200	50 mm	174 mm	200 gr	7,07 oz
ZL-230	50 mm	174 mm	230 gr	8,14 oz