

EO AND STEAM OTOKLAVE TAPE



PRODUCT FEATURES

ZEOSS Ethylene Oxide and Steam Autoclave Sterilization Tape is an adhesive tape used in sterilization processes to show whether a certain gas has been exposed or not. ZEOSS Sterilization Tape has Type 1 indicator.

ZEOSS Ethylene Oxide Sterilization Tapes work by changing color when exposed to certain sterilization conditions.

It is a safe solution for sealing packaging materials.

CODE	PRODUCT NAME	SIZE
ZEOSS EOT1950	ZEOSS Ethylene oxide Sterilization Tape	19mm x 50m
ZEOSS STEAM1950	ZEOSS Steam Sterilization Tape	19mm x 50m