

## HYDROGEN PEROXIDE CHEMICAL INDICATOR



### PRODUCT FEATURES

- Multi-Parameter indicators are designed for two or more of the critical parameters.
- Specifies the excuse and other selected critical parameters during sterilization.
- It is designed for temperature, time and saturated Hydrogen Peroxide control in Hydrogen Peroxide Sterilization process. Thanks to its lamination coating, it prevents possible ink transfer. Thanks to the gradual color change of the indicator, it provides the opportunity to comment on the excuse of the process. The reference color, black, is reached when all critical parameters of steam sterilization have been met. Color change remains stable for at least six months. Designed for inside and outside observation.
- It does not contain heavy metals and is non-toxic. It is marketed in the form of 250 strips/box.
- Conforms to ISO 11140 and EN 867-4 Standards.

CODE	PRODUCT NAME	BOX
ZEO-03	HINROGEN PEROXIDE INDICATOR	250 Pieces