

SEALING MACHINE



PRODUCT FEATURES

Product Code: ZEO 11

The device automatically shuts down.

- The outer and inner surfaces are completely stainless steel, suitable for operation in the operating room environment.
- With its optional removable conveyor, you do not need to hold the packages during the closing process.
- The closing speed can be infinitely adjusted in the range of 0-12 m/s with the pot on it.
- The closure width is 12 mm and it makes a very strong adhesive with barrier in the type of sterilization rolls.
- It can close any size of roll with its rotating wheel system.