

**STEAM BIOLOGICAL INDICATOR**



**PRODUCT FEATURES**

ZEOSS Steam biological indicators have been developed to control the Steam sterilization process. The biological indicator contains bacterial spores embedded in filter paper in a plastic tube. The specially formulated culture medium is contained in a sealed glass ampoule. When the spores germinate, liquid color change in the tube.

ZEOSS Steam indicators give results in 24 hours at 55±2 °C.

CODE	ORGANISM	PERIOD	RESULT	PACKAGE	INCUBATOR
ZB 20	Geobacillus Stearothermophilus	EO	48 hours	100 Pieces	ZDL-100