


ZEISS STERILIZATION REELS



ZEISS medical packaging provides safe and reliable sterile barrier systems for medical devices, pharmaceutical solutions and other healthcare products.

Product Features

- Medical Paper provides maximum safety.
- All imprints are located outside the packing area to prevent ink pigment migration to the product.
- Extra reinforced layered film to avoid tear during opening.
- Superior barrier with 60 gsm or 70 gsm medical kraft paper.
- Non-toxic accurate, process indicators for steam and EO gas conforming to the ISO 11140.
- Formaldehyde (FO) and Hydrogen peroxide (plasma) Indicator available upon request.
- Strong seal-strength, visible sealing lines.
- Clean fibre-free opening.
- Opening direction marked.
- Conformity with EN 868-5.

Product Code	Product Name	Product Size	Package Quantities
MDR050200	ZEISS Medical Sterilization Reel	50mmx200m	16
MDR075200	ZEISS Medical Sterilization Reel	75mmx200m	16
MDR100200	ZEISS Medical Sterilization Reel	100mmx200m	8
MDR150200	ZEISS Medical Sterilization Reel	150mmx200m	4
MDR200200	ZEISS Medical Sterilization Reel	200mmx200m	4
MDR250200	ZEISS Medical Sterilization Reel	250mmx200m	2
MDR300200	ZEISS Medical Sterilization Reel	300mmx200m	2
MDR350200	ZEISS Medical Sterilization Reel	350mmx200m	2
MDR400200	ZEISS Medical Sterilization Reel	400mmx200m	2